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English

Destination and use

The foot switches of series (M)KF serve as actuators for medical devices. They control medical appliances together with the control unit downstreamed by the marketer. The devices are provided with reed contacts, microswitches or HALL sensor and are suitable for applications from 1 mA to approx. 5 A.

Mounting / Wiring

Wire the foot switches according to the specified wire colours/terminal labelling. Contact symbols are shown for a not actuated switch.

Notices

The electrical connection may only be carried out by authorised personnel. The application of a foot-switch bearing the CE mark corresponding to the Medical Devices Act does not exempt the marketer from proving again the conformity of the entire system as per Medical Devices Act (MDD 93/42/EEC)! Despite careful production and inspection as well as redundant design of important components, an entirely safe function at any time cannot be guaranteed (in case of error, a switching insert e. g. cannot switch off properly or a spring may break). The remaining risks have to be taken into consideration concerning the evaluation of the entire system! With use of HALL sensors and reed contacts, magnetic fields can lead to an unwanted influence. Additionally, magnetically conducting materials are to be kept away from the direct foot switch area! In case of application of a shielded control cable, the shield may only be used as functional earth conductor! Subject to technical modifications.

Maintenance / Cleaning

In consideration of the following indications, the steute (M)KF foot switches need only minor maintenance. However, according to ambient conditions and application frequency, we recommend a regular check or maintenance with the following steps:

1. Function test, check actuating elements for proper work.
 2. Check housing and connecting cable for damage and harmful dirt.
- When cleaning the device, use only water and a cloth rag impregnated with a mild cleanser.
 - Do not use cleaning agents which may damage the plastic surfaces, such as detergents, abrasive cleansers or solvent-based cleaners (such as benzine, stain remover).
 - Only clean manually.

Disposal

- Observe national, local and legal regulations concerning disposal.
- Recycle each material separately.

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English

Technical data

Applied standards	EN 60601-1:1990 / IEC 601-1:1988; EN 60601-1:2006 / IEC 60601-1:2005; EN 60601-1-2:2007 / IEC 60601-1-2:2007; EN 60601-2-2:2009 / IEC 60601-2-2:2009; ^{*1} EN 60601-2-22:1996 / IEC 601-2-22:1995; ^{*2} EN 60601-2-43:2010 / IEC 60601-2-43:2010; ^{*1} EN 60529:1991 / IEC 60529:1989
Class	class I as per Medical Devices Act / MDD 93/42/EEC
Enclosure	unbreakable, self extinguishing thermoplast
Switching system	reed contact: 1S, 2S, 1W micro switch: 1PW Hall sensor: HS 0-10 V, 0(4)-20 mA
Degree of protection ^{*3}	IP X1 to IP X8 as per IEC 60529
Termination	firmly attached control line
Switching Voltage	max. 25 VAC / 60 VDC
Current consumption	reed contact: max. 1 A, micro switch: max. 5 A, Hall sensor: max. 60 mA
Switching capacity	reed contact: max. 20 VA, micro switch: max. 1250 VA
Temperature range	storage/transport: -25 °C to +70 °C, operation: -10 °C to +60 °C
Relative humidity	storage/transport: 5 % - 100 %, operation: 10 % - 100 %
Air pressure	storage/transport: 500 hPa - 1100 hPa, operation: 800 hPa - 1060 hPa
Degree of pollution	2 (IP X6 to IP X8) and 3 (IP X1 to IP X5)
Altitude	max. 2000 m

*1 does not apply to foot switch which type of protection is less than IP X7

*2 does not apply to foot switch without protection (bar, protective housing, shroud etc.) against unintentional actuation. The manufacturer of the system has to care for required safety precautions.

*3 Foot switches with protection class IP X1 to IP X5 may not be used in areas, where in NORMAL USE liquids are likely to be found!